

- Is very user friendly
- Works in 145 countries
- Can work indoors and outdoors
- Protects privacy and security of user
- Ensures safety for the patient and gives peace of mind to family members and carers
- Inexpensive to buy and use

**BUY NOW!**  
[www.watchie.com](http://www.watchie.com)



- SOS Panic Button**  
Position data instantly transmitted
- ((G)) Fall Alert**  
Integrated G-Sensor can recognise and trigger a fall-alert
- GeoZones**  
Alert if wearer strays outside defined areas
- i Medical Info Alert**  
Transmits relevant medical info to WATCHIE®
- Privacy**  
Respects wearer's privacy
- A-GPS Localisation & Tracking**  
Advanced CellLocate™ system works indoors and outdoors
- Enhanced Battery Management**  
Innovative Triple-Power-Management® extends standby time
- eSIM Integrated SIM-card**  
Works in 145 countries worldwide
- Comprehensive Reports**  
Gives overview on changing privacy settings and battery level
- Definable Access Rights**  
Different user can be added to the system with different read & write privileges

**LOSTnFOUND AG**  
Soodstrasse 55  
CH-8134 Adliswil  
Switzerland

**International & Switzerland**

T: +41 (0)44 500 40 95  
F: +41 (0)44 500 40 98



**Germany**

T: +49 (0)30 567 947 58  
F: +49 (0)30 567 947 59



**Austria**

T: +43 (0)720 518 700  
F: +43 (0)720 518 701



**UK**

T: +44 2035 149 429  
F: +44 2035 149 430



**LOSTnFOUND (France) SAS**

115 rue Saint Dominique  
F-75014 Paris  
France

**France**

T: +33 (0)1 82 883 997  
F: +33 (0)1 82 883 998



**LOSTnFOUND® (Baltics)**

Kasteheina 4  
EE-51011 Tartu  
Estonia

**Nordics, Baltic States & Russia**

T: +372 666 7291  
F: +372 668 1379



**LOSTnFOUND® (Brasil)**

R. dos Maracanãs 80  
Porto Seguro, BA, CEP 45810-000  
Brasil

**Latin America**

T: +55 (21) 3956 0539  
F: +55 (21) 3956 0541



**Protect your Loved Ones**

**watchie**<sup>®</sup>  
powered by LOSTnFOUND.com

E: [welcome@watchie.com](mailto:welcome@watchie.com)

[www.watchie.com](http://www.watchie.com)



**ADVANCED  
GPS PERSONAL LOCATOR**




[www.watchie.com](http://www.watchie.com)



Home Page



View Map



View Events



View Device Management

WATCHIE<sup>®</sup> is a digital aid for sufferers of dementia. WATCHIE<sup>®</sup> has been developed from the Swiss company, LOSTnFOUND Ltd, which has been very successful in previous years in the business of Vehicles Tracking and Moveable Assets. Now, in conjunction with partners from similar industries, a product has been developed to improve the lives and conditions for sufferers of Alzheimer's and dementia.

## Unique Features of WATCHIE<sup>®</sup>

### Location Positioning Outdoors and Indoors

WATCHIE<sup>®</sup> is the first device in the world to use the CellLocate™ positioning system from u-blox, which intelligently combines the GPS satellite and Mobile phone positioning data together.

### Slimline Design

The design of the WATCHIE<sup>®</sup> is a clever blend of form and function. It incorporates the antenna into the actual wristband, thereby overcoming the notorious Chu-Harrington Effect, which is vital for signal strength.

### Battery Life

This innovative design enables a working battery life of many days, even weeks (depending on set up). Key features of the WATCHIE<sup>®</sup> energy concept are: self activation; efficient management of data transfer; and easily interchanged batteries.

### Privacy

Built in security and privacy functions can be altered to suit the patient and situation.

### Various Alarm Modes

WATCHIE<sup>®</sup> has various modes to raise the alarm. The alarm can be transmitted individually to the various recipients via E-mail or SMS text message:

- On-off switch can be set in a simple manner
- WATCHIE<sup>®</sup> can warn of a fall through it's in-built G-sensor
- With Geo zones WATCHIE<sup>®</sup> can be set up with up to 8 working zones
- Battery power management warns of low battery status
- SOS button
- Altering the activated privacy settings will raise the alarm

### Comprehensive Feed-back

The WATCHIE<sup>®</sup> website gives feed-back on battery status or affected alarms in various ways.

*"For complete self confidence, a person needs freedom and security - this is especially true of dementia patients."*

